

BIODESCONTAMINACIÓN

STERIS BIODECON VHP®

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Julian Pons
Steris Iberia S.A.U.

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Intensive care unit environmental surfaces are contaminated by multidrug-resistant bacteria in biofilms: combined results of conventional culture, pyrosequencing, scanning electron microscopy, and confocal laser microscopy

[H. Hu](#), [K. Johani](#), [J.B. Goebel](#), [A.S.W. Jacombs](#), [A. Almatroudi](#), [G.S. Whiteley](#), [A.K. Dey](#), [S. Jensen](#), [K. Vicker](#) 

Received: March 17, 2015; Accepted: May 26, 2015; Published Online: June 25, 2015

Risk of organism acquisition from prior room occupants: a systematic review and meta-analysis

[B.G. Mitchell](#) , [S.J. Dancer](#), [M. Anderson](#), [E. Dehn](#)

Received: June 8, 2015; Accepted: August 4, 2015; Published Online: August 21, 2015

Publication stage: In Press Accepted Manuscript

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Box 1

Summary of problems associated with conventional cleaning and disinfection

Problems associated with the cleaning/disinfection products include:

- Infectiveness of some agents against some pathogens; for example, many frequently used hospital disinfectants are not effective against *C. difficile* spores and norovirus.^{44,85,134}
- Toxicity to staff or the environment.^{44,46}
- Damage to hospital materials and equipment.⁴⁴
- Susceptibility to interference with organic matter on surfaces.⁸⁵
- Potential for biocide/antibiotic cross-resistance.⁸²

Problems with cleaning/disinfection procedures include:

- Adequate distribution of the active agent, given the challenges of the complex hospital environment.³⁵
- Ensuring correct contact time for the microbial reduction achieved *in vitro*.¹³⁴
- Repeatability of the process depends on the operator.³⁵
- Designation of responsibility for various items, particularly complex portable medical equipment.¹³⁵
- Compliance with protocols/policies from an (often) poorly paid, poorly motivated workforce.¹³⁶
- Inadequate training and education of personnel.¹³⁶
- Inadequate time given to do the job properly.¹³⁶
- Insufficient (or non-existent) cleaning prior to disinfection.⁸⁵
- Incorrect formulation of the disinfectant.^{82,137}
- Contamination of cleaning solutions/materials.^{137,138}
- The effectiveness of conventional cleaning and disinfection is difficult to monitor.¹¹

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From J.A. Otter et al. *Journal of Hospital Infection* 83 (2013) 1-13

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Abbott Laboratories, Abiogen Pharma, Acambis, Acusphere, Alcon, Alkermes, Alpha Therapeutics, Allergan, Altana, Alza, Amersham, Amriya, Aradigm, Asta Medica, Astra-Zeneca, Aventis Pasteur, Aventis Pharma, B.Braun, Bausch & Lomb, Baxter, Bayer, Beaumont Hospital, Dublin, Bigmar Pharma, Biochemie, Biochem Vaccines, Bioglan, Biogen, Biomet, Bioproduct Research, Bioreliance, Biotissue, Boeringher Ingelheim, Bristol-Meyers Squibb, Carbomedics, Cardinal Healthcare, Ciba Vision, Cilag, CBL, Corixa, CP Pharmaceuticals, Depuy CMW, Dey Laboratories, Dr.Kade, Dupont Merck, Eisai, Eli Lilly, Essex, Ethicon, Ferrero, Ferring, Fisons-Aventis, Fort Detrick, Fort Dodge, Fresenius, Fujisawa, GBF, G.C. Hanford, Genetech, Genzyme, Glaxo Sminth Kline, Gruenthal, Gruppo Lepetit, Haupt Pharma, Health Canade, Heinrich Pette Institut, Hemosol, Hikma Pharmaceutical, Hoffmann La Roche, Howmedica, Imperial College London, IVAX Laboratories, Karolinska University, Jackson Laboratories, Janssen, Johnson & Johnson, Knoll, Lancaster Laboratories, Leciva, Lederle laboratories, Life technologies, Lohmann Animal Health, LPI, Mallinkrodt Nuclear laboratories, Marsam Pharmaceutical, Max Plank Institut, Mead-Johnson, Mediva, Medtronic, Megabios, Merck, Merial, Meridian Medical, Monsanto, MBZ, MHH, MSD, MRC, Murinus, Nelson Labs, Nordian, Nordmark, Novartis, Nova Pharmaceutical, Novo Nordisk, Nu-Pharm, Nycomed, OJM Pharmaceutical, Oncotec, Orthivita, Osteotech, OZBT, Patheon, Pasteur-Merieux, Pfizer, Pharmacia & Upjohn, Pierre Fabre, Powderjet, PSI, Purdue Biopharma, Ribbon, Rhone Merieux, Roche Pharmaceutical, Roche diagnostics, RPR-Fisons, RVH Belfast, Sandoz, Sanofi Aventis, Schein Pharmaceutical, Schering-Plough, Sintetica, Solvay Pharma, South Dakota State University, Stowers Institute, SIFI, Sulzer Mitroflow, Systemix, Teva Pharmaceutical, Transgene, Smith & Nephew, Royal College of Surgeons Dublin, Tyco, UKE, Unilever, Union Chimique Belge, University of Bielefeld, University of Cambridge, University of Dresden, University of Durham, University of Erlangen, University of Edinburgh, University of Jena, University of Magdenburg, University of Oxford, University of Tennessee, University of Cork, University of WI, University of Zurich, US Army, Vistakon, Walter Reed Army Institute, Warner Lambert, Wyeth, Zenit.

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Dairy, Food and Environmental Sanitation, Vol. 22, No. 11, Pages 868-873
Copyright© International Association for Food Protection, 6200 Aurora Ave., Suite 200W, Des Moines, IA 50322

Vapor Phase Hydrogen Peroxide Decontamination of Food Contact Surfaces

Gerald McDonnell,* George Grignol, and Kathy Antloga
Research & Development, STERIS Corporation
5960 Heisley Road, Mentor, OH 44060-1834

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2007 Steris named global supplier to Tetra Pak

STERIS and Tetra Pak Sign Joint Development Agreement

Collaborative Efforts Will Leverage Each Company's Core Capabilities

BEBERAGE FILLING STERILIZATION

Tetra Pak tests new sterilisation systems in joint deal

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Control of an outbreak of *Acinetobacter baumannii* infections using vaporized hydrogen peroxide


A. Chmielarczyk^{a,*}, P.G. Higgins^b, J. Wojkowska-Mach^a, E. Synowiec^c,
E. Zander^b, D. Romaniszyn^a, T. Gosiewski^a, H. Seifert^b, P. Heczko^a,
M. Bulanda^a

^a Jagiellonian University Medical College, Krakow, Poland

^b Institute for Medical Microbiology, Immunology and Hygiene, University of Cologne, Cologne, Germany

^c John Paul II Regional Teaching Hospital, Krakow

A. Chmielarczyk et al. / Journal of Hospital Infection 81 (2012) 239–245

APLICACIONES DE STERIS VHP®

BIODESCONTAMINACIÓN

STERIS Defense & Industrial Group Case Study

Anthrax Remediation of U.S. State Department Mail Facility SA-32

The Challenge:

The U.S. State Department Mail Facility in Sterling, Virginia was contaminated with Anthrax spores after accepting mail from the Brentwood Post Office in Washington, D.C. in October 2001. The State Department was forced to abandon and quarantine the facility, and was seeking a safe, environmentally friendly and non-corrosive means of remediation.

The Solution:

STERIS's proprietary vaporized hydrogen peroxide technology (VHP®) was selected as the project decontaminate after the State Department's Technical Working Group (TWG) compared and evaluated available remediation technologies. The group, which included representatives from the State Department, the U.S. Environmental Protection Agency, Centers for Disease Control and Prevention, Occupational Health and Safety Administration, Army Corps of Engineers and the State of Virginia utilized criteria such as efficacy, safety, material compatibility, business continuity, and industry experience with the technology in making their selection.

In use for over a decade with STERIS's pharmaceutical and research customers, the VHP technology was modified by STERIS's Defense & Industrial (D&I) Group for application in large-scale, complex, building interiors. The VHP technology is highly effective against pathogenic organisms and bacterial spores, and unlike many other antimicrobials, it breaks down into water vapor and oxygen. Additionally, the technology will not produce hazardous substances, is environmentally friendly, non-carcinogenic, non-corrosive at use concentration, and does not leave residues.

Because the STERIS VHP system employs a low concentration, dry process, it is compatible with a wide variety of materials, including sensitive electronics, equipment, and a broad range of materials and finishes. It can be used over a wide range of temperatures, works quickly, and because of its scalability, is effective in small or spacious areas.

In addition to providing its VHP technology, the STERIS D&I Group developed a comprehensive remedial action plan and health and safety protocols that were utilized in the successful decontamination of the facility. As

A STERIS Defense & Industrial technician collects biological indicators which were used to validate the efficacy of the decontamination process.

part of its comprehensive initiative, STERIS D&I also drafted and supplied quality and sterility assurance protocols.

The Result:

Remediation of SA-32 occurred between July and August 2003, and was cleared for reuse on December 5, 2003. The approach for decontamination of SA-32 was based on the goal of treating the entire facility section by section, dividing areas into no more than 200,000 cubic feet. These zones were fumigated with VHP as part of a controlled, staged process. In preparation for fumigation, the facility was separated into 10 areas using partitions constructed from polyethylene sheeting with wooden 2' x 4' supporting structures. VHP was pumped into these zones through air intake ducts erected for this purpose.

Measured by real-time VHP concentration monitors, the VHP process met the prescribed specifications and process variables for each of the four phases of the decontamination cycle. Prior to applying the VHP sterilant, chemical and biological indicators were placed throughout the facility. Following VHP application,

chemical indicators exhibited the required colorimetric changes and biological indicators were negative for growth of the indicator organism. Numerous surfaces throughout the facility were swabbed, cultured, and analyzed by methods approved by the Naval Medical Research Center, Silver Springs, Maryland. All 619 samples were negative for *Bacillus anthracis*.

The Environmental Clearance Committee (ECC) stated that the remediation approach used at the facility was thorough in both design and implementation. The ECC achieved its goal of using the best science and technology available to minimize residual risk to the greatest extent possible. The ECC concluded that the decontamination process for the facility as set forth by the TWG and performed by the State Department Project Team, which included STERIS, was safely and effectively completed.

Large scale VHP decontamination system deployed at the U.S. State Department mail processing center in Sterling, Virginia.

To learn more about the STERIS Defense & Industrial Group, visit www.steris.com or call Kevin Marsh at 440-392-7660.

M2712EN 2004-10, Rev. A © 2004 STERIS Corporation. All rights reserved.

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